

TEST REPORT

NASA CR-

141565

FLUID INFUSION SYSTEM

(NASA-CR-141565) FLUID INFUSION SYSTEM
(Beckman Instruments, Inc., Anaheim, Calif.)
15 p HC \$3.25 CSCL 20D


N75-15904

Unclass

G3/34 08946

December 17, 1974

National Aeronautics and Space Administration
Lyndon B. Johnson Space Center
R&T Procurement Branch
Houston, Texas 77058

**Beckman®**

INSTRUMENTS, INC.

ADVANCED TECHNOLOGY OPERATIONS
ANAHEIM, CALIFORNIA 92806

CONTENTS

PARAGRAPH		PAGE
1.0	INTRODUCTION.	1
2.0	ENGINEERING TESTING	1
2.1	Flow Tests	1
2.1.1	Breadboard Flow Tests	1
2.1.2	Prototype Flow Testing.	2
2.2	Bubble Alarm Tests.	7
2.3	Pneumatic System Performance.	7
3.0	ACCEPTANCE TESTING.	7
3.1	Flow Tests	7
3.2	Animal Testing.	9
APPENDIX: ACCEPTANCE TESTING DATA SHEETS		

PRECEDING PAGE BLANK NOT FILMED

1.0 INTRODUCTION

This report describes and summarizes performance testing carried out in the development of the prototype zero-g fluid infusion system under NASA-JSC Contract NAS9-14080. Engineering tests were performed in the course of development, both on the original breadboard device and on the prototype system. This testing was aimed at establishing baseline system performance parameters and facilitating improvements. Acceptance testing was then performed on the prototype system to verify functional performance. Acceptance testing included a demonstration of the fluid infusion system on a laboratory animal.

2.0 ENGINEERING TESTING

The procedures followed and the results obtained in engineering tests are described in sections 2.1 through 2.3.

2.1 Flow Tests

2.1.1 Breadboard Flow Tests

The most important system performance parameter is accuracy of flow regulation. Early in the program, tests were made to establish the flow characteristics produced by the currently most widely used technique--declining hydrostatic head--and the ideal situation of constant fluid pressure applied uniformly over the entire bag exterior. A pressure infusor currently available off the shelf was also tested.

The technique used for measuring flow rate was simply timing the delivery into a graduated cylinder. The volume increment used for this purpose was 5 ml for low flow rates (5 ml/min. nominal), and 50 ml for higher flow rates (50 to 100 ml/min. nominal). Experimental accuracy inherent in this method is on the order of 2 to 3 percent.

Results of the early investigations are shown in Figures 1 and 2. The measured flow rates (Figure 1) and the measured liquid pressure (Figure 2) in the bags are shown as functions of the cumulative liquid volume expressed from the bag. Figure 1 shows two useful baselines. Curve (1) represents a practical ideal--flow is constant except for secondary effects due to elasticity of the fluid bag walls. Curve (2) represents the results obtained in current every-day practice, and shows the substantial decline in flow rate that occurs as the contents of a fluid bag are delivered. Figure 2 shows the flow curves obtained from two devices: the hand-pumped Fenwal pressure infusor, and the original Beckman breadboard fluid infusion system. A factor common to both was the fact that the pneumatic pressure in the device was held constant and it was assumed that this would result in constant liquid pressure. The curves show that this was not the case, however. The Fenwal device did a better job of producing constant flow because of the more conformal nature of the device-to-fluid bag interface. Contact area between bag and pressure plate in the original breadboard was highly variable, and this resulted in a severe decline in flow. Discovery of this effect led to a new conceptual design that was embodied in the prototype fluid infusion system.

2.1.2 Prototype Flow Testing

The superior pneumatic design adopted for the prototype system avoided the gross flow rate errors noted previously. Large flow anomalies were sometimes noted in the prototype engineering tests; however, these were associated with failure modes rather than inherent deficiencies of the design, and they were eliminated. One effect that was found to be systematic, however, was a slow decay in flow at the lower flow rates. This is believed due to elasticity and dimensional stability effects in the flexible tubing. No other consistent effects were noted, and it was found possible to produce a constant flow within $\pm 8\%$ at the lower flow rates and $\pm 5\%$ at the higher flow rates. Some representative flow curves obtained in engineering tests of the prototype system are shown in Figures 3 and 4.

FIGURE 1
FLOW PERFORMANCE--PRESSURE CHAMBER
AND CONVENTIONAL METHOD OF DECREASING HYDROSTATIC HEAD




FIG 1. FLOW RATE VS
VOLUME EXPENDED
(1) BAG IN PRESS. CHAMBER
CONSTANT AT 178 mm Hg
(2) GRAVITY FEED--ELEVATED
TO 33 in. (838 mm)

FIGURE 1




FIGURE 2

FLOW PERFORMANCE--COMMERCIAL PRESSURE INFUSOR
AND INITIAL BECKMAN ATO BREADBOARD DESIGN

FIG. 2--FLOW RATE VS VOLUME EXPENDED

- (1) FENVAL PRESSURE INFUSOR
BLADDER AT CONSTANT PRESSURE
(180 mm Hg)
- (2) BELLOW PRESSURE PLATE AT
CONSTANT PRESSURE (180 mm Hg)


IVI FLOW TESTING


6 DEC 74

RUNS OF 6 DEC 74


- RUN 1 BAG 1 100 cc/M NOMINAL
- △ RUN 1 BAG 2
- ◇ RUN 2 BAG 1 50 cc/M NOMINAL
- ◇ RUN 2 BAG 2
- ▽ RUN 3 BAG 1 50 cc/M NOMINAL
- RUN 3 BAG 2
- RUN 4 BAG 1 5 cc/M NOMINAL

FIGURE 3





IVI Flow Testing - 13 Dec 74



IVI PROTOTYPE SYSTEM - ACCEPTANCE TEST : 16 DEC 74

2.2 Bubble Alarm Tests

Tests were conducted to establish the operating threshold of the bubble/no flow alarm on the modified Carolina Medical Electronics blood flowmeter. Air bubbles were injected into the fluid stream through the septum of the drug injection port using a microsyringe. Bubbles of measured size were injected, and a determination made as to the size bubble that would always trip the alarm, and the size that would trip the alarm about 50% of the time. At the higher flow rates, a larger bubble is required for an alarm due to shorter residence time in the flow probe. Results were as follows:

<u>Nominal Flow Rate</u>	<u>Bubble Size for 0.5 Alarm Probability</u>	<u>Bubble Size for 1.0 Alarm Probability</u>
5 ml/min.	.02 ml	.03 ml
50 ml/min.	.08 ml	0.1 ml
100 ml/min.	.20 ml	0.25 ml

2.3 Pneumatic System Performance

The CO₂-powered pneumatic system was proven out in initial operation of the prototype system. No basic problems of gas bottle icing, materials compatibility, etc., were found. The ability of one 25-gram CO₂ bottle to power the system through a useful number of bag evacuations was verified by test. The lower limit on usable CO₂ bottle pressure was found to be 500 psi--i.e., if bottle pressure is not at or above that pressure, a 2-bag run cannot be undertaken with assurance that there will be enough gas in the system to complete the run. The 25-gram bottle was found to consistently supply enough gas for eight complete runs; i.e., to infuse a total of sixteen 500-ml fluid bags.

3.0 ACCEPTANCE TESTING

Acceptance testing of the unit was performed to demonstrate the performance of the system and its suitability for clinical use.

3.1 Flow Tests

Flow tests were repeated using the engineering flow test procedures. The results obtained are shown in Figure 5, and test data sheets are in the

PRECEDING PAGE BLANK NOT FILMED

Appendix. The flow rate was adjusted upward on the low flow run once steady low flow had been demonstrated.

3.2 Animal Testing

Use of the unit in a clinical environment was demonstrated by infusion of 750 ml of normal saline solution into a laboratory animal. The system was operated normally, with complete delivery from the first bag and automatic switching to the second bag. Flow rates over the entire range were selected during the infusion, with the blood flowmeter being used to monitor flow rate. The test was completed successfully.

PRECEDING PAGE BLANK NOT FILMED

APPENDIX

ACCEPTANCE TESTING DATA SHEETS

DATE: 16 Dec 74

FLOW TEST DATA SHEET - IVE UNIT

NOMINAL FLOW RATE 50 ml/min FLUID 0.9% Saline
BOTTLE PRESSURE: 950 PSI REGULATED PRESSURE: 15 PSI
DELIVERY HEAD: AT TIP 3.5" H₂O AT BAG TOP 3.5" H₂O
DELIVERED VOLUME: AT TRIGGERING 455 ML END BAG 2 445 ML

[illegible]

12/16/77

J. Hammond 16 Dec 74

ORIGINAL PAGE IS
OF POOR QUALITY

LA: 1000019

BOTTLE PRESSURE: 950 PSI REGULATED PRESSURE: 15 PSI

DELIVERED VOLUME: AT TRIGGERING 450 ML END BAG 2 415 ML

12/16/74

John Hammond 16 Dec 74

ORIGINAL PAGE IS
OF POOR QUALITY